EN **EYELID WIPES.**  
**EYELIDS, EYELASHES AND FACE CLEANSER.**

For external use only.

**COMPOSITION**

Terpinen-4-ol, Sodium Hyaluronate, D-Panthenol, Aloe Vera, and purified Water.

**DESCRIPTION**

Demodex mites (Demodex folliculorum and Demodex brevis) are the most common ectoparasites infesting the pilosebaceous unit of the skin. Uncontrolled skin Demodex infestation (demodicosis) has been implicated in several diseases in the skin, including rosacea, with papulopustular skin lesions and perifollicular inflammatory infiltrate. Often, once demodicosis occurs in the face, it may spread and flourish in the eyelids and eyelashes.

**INDICATIONS**

NAVIBLEF® WIPES is a formulation recommended for removal of ocular secretions and debris of eyelids and eyelashes caused by allergic and inflammatory conditions associated with Demodex, Blepharitis, Blefarito-Conjunctivitis, Meibomian gland dysfunction (MGD), Rosacea, burning, itchiness, irritation, redness swelling and other eyelid margin diseases on face. It can also be used for post-cataract surgery. In case of abundant secretions and residues, daily application is recommended.

**PACKAGING**

Box with 20 sterile wipes - [REF] W20NW223

**DIRECTION FOR USE**

**CLOSE YOUR EYES DURING APPLICATION, LET THE SOLUTION DRY UP BEFORE OPENING THE EYES AGAIN (1 mn).**

- Wash, rinse and dry hands and face.
- Tear open the aluminum sachet remove and unfold the wipe.
- Close your eyes tightly and gently wipe the lids using side to side motion for approximately 30 seconds or dab the wipes on the affected area. The solution helps to inhibit the growth of the pathogens and not the scrubbing action. The scrubbing in certain cases may cause patient discomfort.
- Keep the eyes closed for 1 minute to allow to dry the treated surface.
- Do not wash off the dried formulation, let it act for a while. For optimum results do not rinse immediately but wait 5 to 10 minutes before rinsing with running water.
- Once finished cleaning the first eyelid, turn the wipe over and use the opposite side to repeat the procedure with your other eyelid. You can use one wipe for both eyelids unless otherwise directed by your eye specialist.
- Discard wipe after use.

**Recommended NAVIBLEF® WIPES usage according to severity of infestation**

Symptom	Frequency	Duration
Mild - Moderate	1 - 2 times a day	6 - 10 weeks
Moderate - Severe	2 - 3 times a day	10 - 12 weeks

**How to use NAVIBLEF® WIPES**

Natural room temperature for the daily cleansing of face, eyelashes and eyelids.

Warm for thorough cleaning in cases of conjunctivitis, blepharconjunctivitis, sty or occlusion of the lacrimal duct. Immerse the package in lukewarm water for 3 to 5 minutes. Remove the wipe keeping it folded and gently massage the eyelid and eyelashes.

Cold for a decongestant effect in suppurative or inflammation of the eyelids. Place the unopened sachet in the refrigerator for at least 1 hour, remove wipe and hold it onto the area concerned for a few minutes.

**PRECAUTIONS**

- Always wash your hands before and after use.
- Remove contact lenses before the application.
- Keep your eyes closed during the application.
- The product is intended for the application on face, eyelids, and eyelashes.
- Keep out of the sight and reach of children.
- Do not use for children under 3 years.
- Do not use in case of known hypersensitivity to any of the ingredients.
- Some patients may experience a slight stinging sensation, this is normal and should dissipate within a few minutes; if not, rinse with running water or Saline Solution.
- If irritation occurs, stop use immediately. Rinse your eyes with water and contact your eye specialist if the irritation continues.
- Do not use after the expiry date stated on the packaging.
- Do not re-use (risk of contamination and loss of efficacy).
- The used wipes must be disposed in the household waste.

For external use only. Do not apply directly into the eyes.

**STORAGE**

- Store the product in a dry place between 5 and 25°C.

The wipes may have some color variation that does not compromise the quality of the product.

Distributor:  
VISUfarma International,  
2405, route Des Dolines - Drakkar D,  
06560 Valbonne, FRANCE

Patent pending  
NOVAX®PHARMA  
Le Coronado  
20, Av. de Fontvieille  
MC 98000 MONACO  
www.novaxpharma.com

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